



T.b.v. van de deelnemers aan de TOA dag op zaterdag 22 April 2017 in de Schakel te Nijkerk, bieden wij de volledige presentatie aan over voorsegment map interpretatie en sleutelfactoren voor een goed onderzoek.

Daarnaast bieden wij u toegang tot de Ziemer e-learning site. U kunt zich hier gratis registreren, ook als u gebruiker bent van andere apparatuur.

<https://www.e-learning.ziemergroup.com>

De site biedt een verzameling korte video's over zowel de toepassingen (ESCRS lezingen e.d.), maar ook duidelijke uitleg over bijv. het begrip elevatie map en nog veel meer. Van harte aanbevolen.

Meer informatie kunt u krijgen bij Laméris Ootech BV te Ede: 030-6008711 of via ootech@ootech.nl.



Corneal Topography & Tomography

Introduction to measurement parameters and basic map interpretation,
and key factors for obtaining good quality measurements

Gregor Schmid, PhD

Senior Expert Clinical Applications (R&D)
Ziemer Ophthalmic Systems AG, Switzerland



Topography & Tomography



GALILEI



Key Principles



**Dual-Scheimpflug
Tomography**



**Placido Disk
Topography**



Scheimpflug Tomography



Scheimpflug Principle





Single Scheimpflug:

Pros:

- Can measure anterior **and** posterior corneal surface
- Allows 3D view of the cornea
- Works well even in dry eyes
- Permits precise pachymetry **when well centered**
- Provides peripheral data

Cons:


- Does not permit precise pachymetry when decentered
- Precision for curvature is reduced in central area

Dual-Scheimpflug:

Compensation of eye decentration allows precise pachymetry under difficult circumstances (e.g. poor fixation, eye movements, inexperienced operator).



Apparent Corneal Thickness when centered




- When centered, the slit light is perpendicular to the surface
- Apparent thickness with the right and left SF camera are equal



Apparent Corneal Thickness when de-centered

Left SF camera:
Thicker



- The slit light is not perpendicular to the surface
- Apparent thickness is thinner/thicker than at center
- Averaging automatically corrects de-centration




Placido Topography





Placido Rings

SF Cameras




Placido Topography



Placido based topographers work on the principle of assessing the reflection of a concentric set of white rings from the convex anterior surface of the cornea



Placido Topography



Placido Topography




Pros:

- Allows accurate measurement of the anterior surface curvature

Cons:

- Suffers from dependency on tear film quality
- Provides no information on posterior surface and pachymetry
- Is limited to central area



Dual Scheimpflug + Placido


- **Combines the best of both worlds to produce the most complete data set of the anterior segment**




Maps and measurements aligned to apex (1st Purkinje images)

Motion Compensation using iris pattern

- a) Lateral motion correction (x/y-directions)
- b) Rotational correction (around z-axis)



a)



b)




HOW TO GET GOOD MEASUREMENTS




1. Alignment

Left-Right/Up-Down




In-Out






Alignment: Purkinje Image 1

Alignment to P1 ↔ Visual Axis




Alignment to Apex ↔ ??






2. Tear Film






3. Partial Lid Closure





4. Eye Movement





Tips for good quality measurements

Make sure...

- Patient is well instructed
- Head is correctly aligned (fore-front, chin, head-turn)
- Mouth/jaw is closed
- Eye is wide open during scans

- Procedure:**
- 1. Align the device**
 - 2. Tell patient to blink one more time and open wide**
 - 3. Fine-align and initiate the measurement quickly**
 - 4. Start again if needed, there is no need to rush**



Measurement Quality Values

	Actual	Recomm.
Motion Comp. ✓	100.0%	85.0%
Placido ✓	93.0%	85.0%
Scheimpflug ✓	98.7%	90.0%
Motion Distance ✓	100.0%	70.0%
Overall Quality	97.5%	

Quick Guide to taking successful Measurements



1. Position Patient

Adjust table height to maximize patient's comfort. Patient's head should be straight, positioned centrally with forehead against the headrest and chin touching the front edge of the chinrest.




2. Adjust Headrest

Align patient's lower lid to black marker on the headrest post. Turn the headrest to the left to measure OD and to the right to measure OS.



3. Dim Room Lights

Dim lighting in the room to prevent interfering light reflections during the capture.



4. Gross Alignment

Enter scan mode.

- Move the device to the OS or OD sign on the table (consider the correct position of your headrest, see 2.)
- Using the joystick, align the measurement head so that the horizontal blue light from the slit lamp is centered on the patient's pupil.
- Move the joystick right or left for horizontal and rotate up or down for vertical alignment. Move forward or backward to put the eye into focus.



5. Fixate Target

Ask patient to fixate on visible red bull's eye target.



6. Center

Align the red crosshair to the 4 white dots in the top view image by rotating the joystick.



7. Fine Alignment

Align the red reference line to the anterior surface of the cornea by moving the joystick forward or backward. The red line should not intersect the cornea, but just touch it.



8. Blink and Open Wide

Ask the patient to blink and keep the eye open wide.



9. Scan



Press the scan button on the joystick.

Quality Check Guide





1 Check quality of Scheimpflug/ Placido ring images

- A** Good quality measurement: measurement is centered and cornea aligned to the red reference line
- B** Poor quality measurement: nose shadow, here at the bottom right, is covering the eye
- C** Poor quality measurement: eye lashes or lid are covering the eye



2

- A** Good quality measurement: cornea is touching the red reference line
- B** Poor quality measurement: cornea is intersecting the red reference line
- C** Poor quality measurement: cornea is too far from the red reference line




3

- A** Head should be straight, positioned centrally with forehead against the headrest and chin touching the front edge of the chinrest. Patient is focusing well on the red reference target.
- B** Head is not positioned correctly; head is tilted, forehead is not touching the headrest, chin is not touching the front edge of the chinrest. Patient is not focusing well on the red reference target. Adjust height of table and chinrest to correct your patient's position.









Anterior Axial Curvature





WTR Astigmatism

- Database
- Measure
- Analyze
- Verify
- Refractive
- Asymmetries
- CLMIX
- Thickness Pro
- PTA
- Wavefront
- IOL Power
- Map x 1
- Mapx4 1
- Mapx4 2
- IOL Calculator
- ZLASIK.
- Comparison
- Densitometry
- Eye Metrics



SimK			
SimK	46.95 D	R	7.19 mm
Flat SimK	45.85 D	174°	R1 7.36 mm
Steep SimK	48.05 D	84°	R2 7.02 mm
Astig	2.20 D	84°	e ² (-Q) 0.34

Anterior Axial Curvature Zones			
Central	47.03 D	7.18 mm	
Mid	46.19 D	7.31 mm	
Periph	43.55 D	7.75 mm	
Kmax	48.68 D	6.93 mm	location x,y 0.14 mm 0.89 mm

Posterior Axial Curvature			
Mean K	-6.70 D	R	5.97 mm
Flat K	-6.57 D	1°	R1 6.09 mm
Steep K	-6.83 D	91°	R2 5.85 mm
Astig	-0.27 D	91°	e ² (-Q) 0.68

Pachymetry			
o Thinnest	620 µm	x,y	0.23 mm -0.45 mm
Central	630 µm	CCT	621 µm
Mid	667 µm		
Periph	725 µm	Corneal Vol.	34.8 mm ³

Total Corneal Power IOL (Ray Traced)			
Mean TCPIOL	45.97 D	Central	45.89 D
Flat TCPIOL	44.84 D	172°	Mid 46.60 D
Steep TCPIOL	47.10 D	82°	Periph 45.66 D
Astig	2.26 D	82°	





Anterior Chamber and Biometry			
WTW, N-T	11.52 mm	Mean Angle	30.9 °
ACV	90 mm ³	Kappa Dist	0.34 mm
ACD	3.05 mm	ASL endo	n/a
+ Pupil Diam	2.71 mm	location x,y	0.29 mm -0.18 mm

Corneal Shape Asymmetry			
KPI	9.4 %	Kprob	8.1 %
CLMIIa	0.94 D	PPK	1.2 %



ATR Astigmatism

- Database
- Measure
- Analyze
- Verify
- Refractive
- Keratoconus
- Decision Tree
- CLMIX
- PTA
- Wavefront
- IOL Power
- Map x 1
- IOL Calculator
- Z LASIK.
- Comparison
- Densitometry
- Eye Metrics







SimK			
SimK	40.65 D	R	8.30 mm
Flat SimK	40.18 D	68°	R1 8.40 mm
Steep SimK	41.12 D	158°	R2 8.21 mm
Astig	0.94 D	158°	e ² (-Q) 0.20
Anterior Axial Curvature Zones			
Central	40.61 D	8.31 mm	
Mid	40.25 D	8.39 mm	
Periph	39.23 D	8.60 mm	
Kmax	41.21 D	8.19 mm	location x,y -0.96 mm 0.53 mm
Posterior Axial Curvature			
Mean K	-5.67 D	R	7.05 mm
Flat K	-5.61 D	161°	R1 7.13 mm
Steep K	-5.74 D	71°	R2 6.97 mm
Astig	-0.13 D	71°	e ² (-Q) 0.58
Pachymetry			
o Thinnest	540 μ m	x,y	0.76 mm -0.79 mm
Central	552 μ m	CCT	545 μ m
Mid	584 μ m		
Periph	622 μ m	Corneal Vol.	30.4 mm ³
Total Corneal Power IOL (Ray Traced)			
Mean TCPIOL	39.84 D	Central	39.76 D
Flat TCPIOL	39.24 D	Mid	40.43 D
Steep TCPIOL	40.45 D	Periph	41.06 D
Astig	1.21 D	159°	
Anterior Chamber and Biometry			
WTW, N-T	12.52 mm	Mean Angle	25.7 °
ACV	100 mm ³	Kappa Dist	0.39 mm
ACD	2.79 mm	ASL endo	n/a
+ Pupil Diam	3.31 mm	location x,y	0.16 mm -0.35 mm
Keratoconus Probability			
KPI	0.0 %	Kprob	0.3 %
CLMIIa	0.08 D	PPK	0.2 %



Oblique Astigmatism

- Database
- Measure
- Analyze
- Verify
- Refractive
- Keratoconus
- Decision Tree
- CLMIX
- PTA
- Wavefront
- IOL Power
- Map x 1
- IOL Calculator
- ZLASIK.
- Comparison
- Densitometry
- Eye Metrics



Simk			
Simk	41.23 D		R 8.19 mm
Flat Simk	41.14 D	146°	R1 8.20 mm
Steep Simk	41.33 D	56°	R2 8.17 mm
Astig	0.20 D	56°	e ² (-Q) 0.16

Anterior Axial Curvature Zones			
Central	41.24 D	8.18 mm	
Mid	40.85 D	8.26 mm	
Periph	40.21 D	8.39 mm	
Kmax	41.85 D	8.06 mm	location x,y 0.27 mm 1.27 mm

Posterior Axial Curvature			
Mean K	-5.88 D		R 6.80 mm
Flat K	-5.73 D	8°	R1 6.97 mm
Steep K	-6.02 D	98°	R2 6.64 mm
Astig	-0.29 D	98°	e ² (-Q) 0.26

Pachymetry			
o Thinnest	527 μm		x,y -0.19 mm -0.06 mm
Central	535 μm		CCT 527 μm
Mid	584 μm		
Periph	654 μm		Corneal Vol. 30.6 mm ³


Total Corneal Power IOL (Ray Traced)			
Mean TCPIOL	40.29 D		Central 40.26 D
Flat TCPIOL	40.11 D	116°	Mid 40.69 D
Steep TCPIOL	40.46 D	26°	Periph 41.93 D
Astig	0.35 D	26°	

Anterior Chamber and Biometry			
WTW, N-T	12.43 mm		Mean Angle n/a
ACV	n/a		Kappa Dist 0.13 mm
ACD	n/a		ASL endo n/a
+ Pupil Diam	2.63 mm		location x,y 0.03 mm 0.13 mm

Keratoconus Probability			
KPI	0.0 %		Kprob 1.1 %
CLMIIaa	0.90 D		PPK 1.1 %




Curvature



$$\text{Curvature} = n/r$$

Anterior:

- SimK = calculated with keratometric index (1.3375) assuming constant ratio of anterior & posterior cornea




- Scale with adjustable step-size
- Towards red/white: steeper, towards blue = flatter

Posterior:

- Posterior K = calculated with index of refraction of cornea (1.376) and aqueous humour (1.336) with ray-tracing
- Negative power



Curvature




- **Instantaneous Curvature** = local representation of curvature (C_1, C_2)
- **Axial Curvature** = smoothened curvature (radius extended to reference axis: A_1, A_2), making steep areas flatter, and flat areas steeper



Instantaneous Curvature


- Database
- Measure
- Analyze
- Verify
- Refractive
- Asymmetries
- CLMIX
- Thickness Pro
- PTA
- Wavefront
- IOL Power
- Map x 1
- Mapx4 1
- Mapx4 2
- IOL Calculator
- Z LASIK.
- Comparison
- Densitometry
- Eye Metrics



SimK			
SimK	46.95 D	R	7.19 mm
Flat SimK	45.85 D	174°	R1 7.36 mm
Steep SimK	48.05 D	84°	R2 7.02 mm
Astig	2.20 D	84°	e ² (-Q) 0.34
Posterior Axial Curvature			
Mean K	-6.70 D	R	5.97 mm
Flat K	-6.57 D	1°	R1 6.09 mm
Steep K	-6.83 D	91°	R2 5.85 mm
Astig	-0.27 D	91°	e ² (-Q) 0.68
Total Corneal Power IOL (Ray Traced)			
Mean TCPIOL	45.97 D	Central	45.89 D
Flat TCPIOL	44.84 D	172°	Mid 46.60 D
Steep TCPIOL	47.10 D	82°	Periph 45.66 D
Astig	2.26 D	82°	
Pachymetry			
o Thinnest	620 μm	x,y	0.23 mm -0.45 mm
Central	630 μm	CCT	621 μm
Mid	667 μm		
Periph	725 μm	Corneal Vol.	34.8 mm ³
Anterior Instantaneous Curvature Zones			
Central	46.80 D	7.21 mm	
Mid	43.21 D	7.81 mm	
Periph	33.90 D	9.96 mm	
Kmax	48.85 D	6.91 mm	location x,y 0.15 mm 0.68 mm
Anterior Chamber and Biometry			
WTW, N-T	11.52 mm	Mean Angle	30.9 °
ACV	90 mm ³	Kappa Dist	0.34 mm
ACD	3.05 mm	ASL endo	n/a
+ Pupil Diam	2.71 mm	location x,y	0.29 mm -0.18 mm
Corneal Shape Asymmetry			
KPI	9.4 %	Kprob	8.1 %
CLMIIa	0.94 D	PPK	1.2 %



Instantaneous Curvature: mm




SimK			
SimK	43.94 D	R	7.68 mm
Flat SimK	43.70 D	53°	R1 7.72 mm
Steep SimK	44.18 D	143°	R2 7.64 mm
Astig	0.48 D	143°	e ² (-Q) 0.17
Posterior Axial Curvature			
Mean K	-5.90 D	R	6.78 mm
Flat K	-5.78 D	110°	R1 6.92 mm
Steep K	-6.01 D	20°	R2 6.65 mm
Astig	-0.23 D	20°	e ² (-Q) -1.36
Total Corneal Power IOL (Ray Traced)			
Mean TCPIOL	43.28 D	Central	43.30 D
Flat TCPIOL	42.95 D	45°	Mid 43.98 D
Steep TCPIOL	43.61 D	135°	Periph 45.64 D
Astig	0.67 D	135°	
Pachymetry			
o Thinnest	450 μm	x,y	0.66 mm -1.35 mm
Central	459 μm	CCT	455 μm
Mid	489 μm		
Periph	549 μm	Corneal Vol.	25.6 mm ³
Anterior Instantaneous Curvature Zones			
Central	43.76 D		
Mid	43.25 D		
Periph	41.70 D		
Kmax	44.43 D	7.60 mm	location x,y -0.08 mm 0.06 mm
Anterior Chamber and Biometry			
WTW, N-T	11.23 mm	Mean Angle	52.5 °
ACV	233 mm ³	Kappa Dist	0.07 mm
AQD	3.93 mm	ASL endo	n/a
+ Pupil Diam	6.22 mm	location x,y	0.03 mm 0.06 mm
Corneal Shape Asymmetry			
CLMIIa	0.10 D	PPK	0.2 %

Show Eye Transparency: Map Eye

Instantaneous Curvature: mm

Posterior Instantaneous Curvature Radius [mm]

No map alignment



SimK n.1.3375


	SimK 43.94 D		R	7.68 mm
	Flat SimK 43.70 D	53°	R1	7.72 mm
	Steep SimK 44.18 D	143°	R2	7.64 mm
	Astig 0.48 D	143°	e ² (-Q)	0.17
	Posterior Axial Curvature			
	Mean K -5.90 D		R	6.78 mm
	Flat K -5.78 D	110°	R1	6.92 mm
	Steep K -6.01 D	20°	R2	6.65 mm
	Astig -0.23 D	20°	e ² (-Q)	-1.36
	Total Corneal Power IOL (Ray Traced)			
	Mean TCPIOL 43.28 D		Central	43.30 D
	Flat TCPIOL 42.95 D	45°	Mid	43.98 D
	Steep TCPIOL 43.61 D	135°	Periph	45.64 D
	Astig 0.67 D	135°		
	Pachymetry			
	o Thinnest	450 μm	x,y	0.66 mm -1.35 mm
	Central	459 μm	CCT	455 μm
	Mid	489 μm		
	Periph	549 μm	Corneal Vol.	25.6 mm ³
	Posterior Instantaneous Curvature Zones			
	Central	-6.03 D		
	Mid	-5.93 D		
	Periph	-5.73 D		
	Anterior Chamber and Biometry			
	WTW, N-T	11.23 mm	Mean Angle	52.5 °
	ACV	233 mm ³	Kappa Dist	0.07 mm
	AQD	3.93 mm	ASL endo	n/a
	+ Pupil Diam	6.22 mm	location x,y	0.03 mm 0.06 mm
	Corneal Shape Asymmetry			
	CLMIaa	0.10 D	PPK	0.2 %

Show Eye Transparency: Map

Eye




Refractive Power




- Calculated by ray-tracing through the anterior corneal surface

$$\text{Power} = n/f$$

- f = focal length, n = refractive index
- Focal length is determined as the distance from the reference plane to the intersection of the ray with the central axis
- To determine the focal length the reference plane is the anterior corneal surface in this case



Total Corneal Power (TCP)




$$TCP = n/f$$

- Calculated by ray-tracing through the anterior and posterior corneal surface, taking into account the actual indices of refraction (n)


- TCP_{IOL} is calculated using $n_{aqueous} = 1.336$. To determine the focal length (f), the reference plane is the posterior corneal surface
- TCP_{IOL} gives a more realistic characterisation of the anterior and posterior corneal surface.

The results will however deviate from SimK and cannot be used in traditional IOL formulas.







Pachymetry



- Shows corneal thickness profiles in 20 μm steps
- Towards red/white: thinning, towards blue: thickening
- Thinnest point = indicated by a small circle
- CCT = central corneal thickness; corresponds to central value of the map



Elevation






- Elevation requires a reference (plane, sphere, asphere,...). **GALILEI: BFS, BFA, BFTA**
- Example: Elevation of a mountain





Elevation







- : Cornea
- - - : Best-Fit-Sphere (BFS)



Corneal Shape Asymmetries


- Database
- Measure
- Analyze
- Verify
- Refractive
- Asymmetries
- CLMIX
- Thickness Proj
- PTA
- Wavefront
- IOL Power
- Map x 1
- Mapx4 1
- Mapx4 2
- IOL Calculator
- ZLASIK.
- Comparison
- Densitometry
- Eye Metrics



SimK				n 1.3375	
SimK	44.73 D		R	7.55 mm	
Flat SimK	44.41 D	158°	R1	7.60 mm	
Steep SimK	45.05 D	68°	R2	7.49 mm	
Astig	0.64 D	68°	e ² (-Q)	0.53	
Anterior Instantaneous Curvature (Tangential)				n 1.3375	
Mean K	44.13 D		R	7.65 mm	
Flat K	43.83 D	159°	R1	7.70 mm	
Steep K	44.42 D	69°	R2	7.60 mm	
Astig	0.59 D	69°	e ² (-Q)	0.53	
Total Corneal Power IOL (Ray Traced)					
Mean TCPIOL	43.86 D		Central	43.75 D	
Flat TCPIOL	43.59 D	151°	Mid	43.66 D	
Steep TCPIOL	44.13 D	61°	Periph	42.34 D	
Astig	0.54 D	61°			
Anterior Chamber and Biometry					
WTW, N-T	12.42 mm		Mean Angle	24.8 °	
ACV	81 mm ³		Kappa Dist	0.25 mm	
ACD	2.85 mm		ASL endo	n/a	
+ Pupil Diam	1.90 mm		location x,y	-0.15 mm 0.20 mm	
Pachymetry					
o Thinnest	608 µm		x,y	-0.34 mm -0.49 mm	
Central	618 µm		CCT	611 µm	
Mid	656 µm				
Periph	712 µm		Corneal Vol.	34.3 mm ³	
Corneal Shape Asymmetry & Indices					
KPI	28.1 %		Kprob	58.9 %	
CLMIIa	1.33 D		PPK	2.6 %	
I-S	0.83 D		DSI	2.48 D	
SAI	1.39 D		OSI	2.31 D	
SRI	1.02 D		CSI	0.88 D	
IAI	0.50 D		ACP	44.83 D	
AA	100.0 %		SDP	2.06 D	









Total Corneal Wavefront



Corneal Wavefront Aberrations:

- Path length differences between the actual wavefront and a plane wavefront at the entrance pupil, normally expressed in μm
- Most common aberrations: spherical aberration, astigmatism, coma, defocus



- **Spherical Aberration:** occurs when light experience stronger refractive power at the periphery of the cornea, resulting in a region of defocused light and decreased image quality.

Comparison





GALILEI File Examination Options

Schmid Gregor
2016 Mar 24

A OS/OD OD 2016 Aug 23 - 08:32 Beta-Messung
B OS/OD OD 2016 Aug 18 - 11:52

DAN EMR  

- Database
- Measure
- Analyze
- Verify
- Refractive
- Keratoconus
- Decision Tree
- CLMIX
- PTA
- Wavefront
- IOL Power
- Map x 1
- IOL Calculator
- ZLASIK
- Comparison
- Densitometry
- Eye Metrics
- Settings



A SimK n 1.3375

SimK	44.19 D		R	7.64 mm
Flat SimK	43.59 D	176°	R1	7.74 mm
Steep SimK	44.79 D	86°	R2	7.54 mm
Astig	1.20 D	86°	e ² (-Q)	0.16

Anterior Chamber and Biometry

WTW, N-T	12.34 mm	Mean Angle	n/a
ACV	143 mm ³	Kappa Dist	0.03 mm
AQD	3.55 mm	ASL endo	n/a
+ Pupil Diam	3.81 mm	location x,y	0.00 mm 0.03 mm

B SimK n 1.3375

SimK	42.22 D		R	7.99 mm
Flat SimK	41.88 D	166°	R1	8.06 mm
Steep SimK	42.56 D	76°	R2	7.93 mm
Astig	0.67 D	76°	e ² (-Q)	0.23

Anterior Chamber and Biometry

WTW, N-T	11.91 mm	Mean Angle	30.7°
ACV	92 mm ³	Kappa Dist	0.19 mm
AQD	2.75 mm	ASL endo	n/a
+ Pupil Diam	4.57 mm	location x,y	-0.13 mm 0.14 mm

Difference A-B merged_bfs_elevation_acs Zone Differences


	Mean	RMS
Central	0.00 D	0.00 D
Mid	-0.00 D	0.00 D
Periph	0.01 D	0.02 D
Total	n/a	0.01 D



Eye Metrics



- Database
- Measure
- Analyze
- Verify
- Refractive
- Keratoconus
- Decision Tree
- CLMIX
- PTA
- Wavefront
- IOL Power
- Map x 1
- IOL Calculator
- Z LASIK.
- Comparison
- Densitometry
- Eye Metrics

Viewing Angle



OD OS

1 / 17



Filters

Inversed Image No enhancement

Colorized Image Cornea focussed enhancement

SimK				n 1.3375
SimK	43.15 D		R	7.82 mm
Flat SimK	42.75 D	168°	R1	7.89 mm
Steep SimK	43.56 D	78°	R2	7.75 mm
Astig	0.81 D	78°	e ² (-Q)	-0.12
Posterior Axial Curvature				
Mean K	-6.35 D		R	6.29 mm
Flat K	-6.11 D	5°	R1	6.54 mm
Steep K	-6.60 D	95°	R2	6.06 mm
Astig	-0.48 D	95°	e ² (-Q)	0.68
Anterior Axial Curvature Zones				n 1.3375
Central	43.04 D			7.84 mm
Mid	43.24 D			7.81 mm
Periph	43.17 D			7.82 mm
Kmax	44.66 D		location x,y	0.00 mm 2.90 mm
Total Corneal Power IOL (Ray Traced)				
Mean TCPIOL	42.03 D		Central	41.81 D
Flat TCPIOL	41.74 D	152°	Mid	43.47 D
Steep TCPIOL	42.32 D	62°	Periph	46.28 D
Astig	0.58 D	62°		
Pachymetry				
o Thinnest	536 μm		x,y	-0.97 mm -0.52 mm
Central	556 μm		CCT	544 μm
Mid	606 μm			
Periph	633 μm		Corneal Vol.	31.3 mm ³
Anterior Chamber and Biometry				
WTW, N-T	11.83 mm		Mean Angle	30.4 °
ACV	104 mm ³		Kappa Dist	0.40 mm
ACD	3.00 mm		ASL endo	n/a
+ Pupil Diam	2.70 mm		location x,y	-0.40 mm -0.05 mm
Keratoconus Probability				
KPI	0.0 %		Kprob	1.9 %
CLMIaa	0.92 D		PPK	1.1 %



Densitometry

Database

Measure

Analyze

Verify

Refractive

Keratoconus

Decision Tree

CLMIX

PTA

Wavefront

IOL Power

Map x 1

IOL Calculator


ZLASIK.

Comparison

Densitometry


Eye Metrics

Viewing Angle



OD OS

1 / 17



- Filters
- Inversed Image
 - No enhancement
 - Colorized Image
 - Cornea focussed enhancement

SimK				n 1.3375
SimK	43.15 D		R	7.82 mm
Flat SimK	42.75 D	168°	R1	7.89 mm
Steep SimK	43.56 D	78°	R2	7.75 mm
Astig	0.81 D	78°	e ² (-Q)	-0.12
Posterior Axial Curvature				
Mean K	-6.35 D		R	6.29 mm
Flat K	-6.11 D	5°	R1	6.54 mm
Steep K	-6.60 D	95°	R2	6.06 mm
Astig	-0.48 D	95°	e ² (-Q)	0.68
Anterior Axial Curvature Zones				n 1.3375
Central	43.04 D			7.84 mm
Mid	43.24 D			7.81 mm
Periph	43.17 D			7.82 mm
Kmax	44.66 D		location x,y	0.00 mm 2.90 mm
Total Corneal Power IOL (Ray Traced)				
Mean TCPIOL	42.03 D		Central	41.81 D
Flat TCPIOL	41.74 D	152°	Mid	43.47 D
Steep TCPIOL	42.32 D	62°	Periph	46.28 D
Astig	0.58 D	62°		
Pachymetry				
o Thinnest	536 μm		x,y	-0.97 mm -0.52 mm
Central	556 μm		CCT	544 μm
Mid	606 μm			
Periph	633 μm		Corneal Vol.	31.3 mm ³
Anterior Chamber and Biometry				
WTW, N-T	11.83 mm		Mean Angle	30.4 °
ACV	104 mm ³		Kappa Dist	0.40 mm
ACD	3.00 mm		ASL endo	n/a
+ Pupil Diam	2.70 mm		location x,y	-0.40 mm -0.05 mm
Keratoconus Probability				
KPI	0.0 %		Kprob	1.9 %
CLMIaa	0.92 D		PPK	1.1 %

Clinical Benefits GALILEI



- Reliable and fast topography and tomography screening
- Highly accurate anterior and posterior curvature assessment for sensitive keratoconus screening
- Maps and data aligned to the same reference – the visual axis
- Spherical and aspherical aberrations for wave front guided treatments and toric IOL selection



GALILEI G4

Dual Scheimpflug Analyzer

GALILEI G6

Lens Professional